







Date: N°: Patient: Age: <input type="checkbox"/> Male <input type="checkbox"/> Female Weight: Height: Doctor: Centre: 1 st Order: <input type="checkbox"/> Replacement: <input type="checkbox"/>	Vest: N° of units: <input type="checkbox"/>  Without sleeves: 9101 <input type="checkbox"/> Short sleeves: 9111 <input type="checkbox"/> Long sleeves: 9121 <input type="checkbox"/>	Entire trunk: N° of units: <input type="checkbox"/>  Without sleeves: 9101 <input type="checkbox"/> Short sleeves: 9111 <input type="checkbox"/> Long sleeves: 9121 <input type="checkbox"/>	Arm sleeve: N° of units: <input type="checkbox"/>  Ref. 9201 R <input type="checkbox"/> L <input type="checkbox"/>	Arm sleeve and shoulder Flap N° of units: <input type="checkbox"/>  Ref. 9202 R <input type="checkbox"/> L <input type="checkbox"/>
---	--	--	---	--



BASIC MEASUREMENTS		R	L
1	Circumference of neck		
2	Circumference of chest under armpit		
3	Circumference of under breast perimeter		
4	Circumference of waist		
5	Circumference of hips		
6	Circumference of arm		
7	Circumference of elbow		
8	Circumference of 1/2 forearm		
9	Circumference of wrist		
10	Armpit perimeter		
11	Distance from shoulder to waist (front)		
12	Distance from shoulder to breast		
13	Inside length to arm		
14	Outside length to arm		
15	Length from wrist to elbow		
16	Width of shoulder		
17	Distance from shoulder to waist (rear)		
18	Distance from shoulder to end of garment		

ADDITIONAL MEASUREMENTS		
LENGTH (from wrist) each 5 cms.	ARM CIRCUMFERENCE	
	R	L

R = right
L = left

- Fastening Front
 Rear
 Velcro
 Zip
- Very open
 High neck Height
- ADDITIONAL MEASUREMENT FOR WOMEN
 Perimeter on nipple line
- Brassière size

COMMENTS